

CytoFLEX LX

Flow Cytometer

UV-Violet-Blue-Yellow Green- Red-Infrared

(U-V-B-Y-R-I) Series



OPTICS

EXCITATION OPTICS

The instrument has the capacity for 23 parameters, including 21 for fluorescence detection. The fully activated instrument includes three fluorescent channels from the 355 nm (UV) laser, five from the 405 nm (Violet) laser, three from the 488 nm (Blue) laser, five from the 561 nm (Yellow Green) laser, three from the 638 nm (Red) laser, and two from the 808 nm (Infrared) laser. Instruments with as few as 14 fluorescent channels activated are available with the ability to activate additional parameters as needed by purchasing an activation key.

LASER SPECIFICATIONS

Spatially Separated Laser Options (Beam Spot Size: 5 µm x 80 µm)

Laser	Wavelength	Power
UV	355 nm	20 mW
Violet	405 nm	80 mW
Blue	488 nm	50 mW
Yellow Green	561 nm	30 mW
Red	638 nm	50 mW
Infrared	808 nm	60 mW

FLOW CELL

Fixed integrated optics and quartz flow cell design with >1.3 numerical aperture.

Flow Cell dimensions: 420 µm x 180 µm internal diameter

FORWARD SCATTER DETECTION

Proprietary Homodyne FSC sensor system using silicon photodiodes with built in 488/8 nm bandpass filter.

BANDPASS FILTERS

Includes 23 repositionable filters

405/10	675/30 (2)
405/30	690/50
450/45 (2)	710/50
525/40 (3)	712/25
585/42	763/43 (3)
610/20 (3)	840/20
660/10 (2)	885/40

FLUORESCENCE AND SIDE SCATTER DETECTION

Fluorescence and side scatter light delivered by fiber optics to Avalanche Photo Diode detector arrays. Proprietary design ensures high performance, high efficiency, low-noise signal detection. Emission profiles are collected using reflective optics and single transmission band pass filters.

VIOLET SIDE SCATTER CONFIGURATION

Option to configure Avalanche Photo Diode detector array to collect side scatter signal from Violet (405 nm) laser. The configured channel (VSSC) can be used to better resolve nanoparticles.

QUALITY CONTROL

For detection channel off of the 405, 488, 561, and 638 nm laser, CytExpert QC automation pass/fail criteria is rCV ≤5.0%. For detection channels off of the 355 nm and 808 nm laser, the criteria is ≤7.0%.

PERFORMANCE

SCATTER RESOLUTION

Blue (488 nm) Side Scatter Resolution: <300 nm

Violet (405 nm) Side Scatter Resolution (VSSC):
80 nm polystyrene particles

Scatter performance is optimized for resolving human lymphocytes, monocytes, and granulocytes as well as nanoparticles.

CARRYOVER

Single Tube Format: ≤1.0%

Plate Loader Format: <0.5%

SENSITIVITY

FITC: <30 molecules of equivalent soluble fluorochrome (MESF-FITC) from the 488 nm laser.

PE: <10 molecules of equivalent soluble fluorochrome (MESF-PE) from the 488 nm laser.

FLUORESCENCE RESOLUTION

The CytoFLEX Flow Cytometer is capable of achieving <3% rCV with alignment verification particles capable of rCVs <3%.

ELECTRONICS

NOMINAL ACQUISITION RATE

30,000 events per second with all configured parameters

Software capability to modify window extension parameter and to control abort rate during high event rate signal processing

SIGNAL PROCESSING

Fully digital system with 7 decade data display

SIGNAL

Pulse area, height for every channel, width for one selectable channel

FLUIDICS

Ultra-low pressure peristaltic sheath and sample delivery system

Sheath Fluid Filter and Sample Pump Tubing can be replaced by the user (no service visit required)

SAMPLE FLOW RATES

Fixed Flow Rates: 10, 30 and 60 $\mu\text{L}/\text{min}$

Custom Flow Rate Control mode from 10 to 240 $\mu\text{L}/\text{min}$ in 1 μL increments.

FLUID CAPACITY

10L cubitainers

AUTOMATED MAINTENANCE FUNCTIONS

System Startup, Sample Mixing, Backflush, Prime, Shutdown, Deep Clean

SAMPLE INPUT FORMATS

5 mL (12 x 75 mm) polystyrene and polypropylene tubes

1.5 mL and 2 mL microcentrifuge tubes

PLATE LOADER FORMATS

96-well Standard Flat, U and V bottom plates, and 96-deep well plates. Refer to CytoFLEX Plate Loaders Technical Specification Sheet for details on all plate loader options.

DATA MANAGEMENT

SOFTWARE

The CytExpert software is a full-feature software package that controls the instrument's operation, collection of experiment data, and analysis of the results.

Three different installation modes are available, depending on the level of security required.

The Default installation requires no user login.

For multiuser instruments, the User Management mode requires user login and contains features for role management.

Electronic Records Management installation provides tools that facilitate compliance with 21 CFR Part 11, Electronic Records and Electronic Signatures.

Gravimetric calibration for absolute counts within CytExpert Software.

An API (Application Programming Interface) is available and allows external software to perform operations such as running methods and for basic control of the plate loader.

If desired, export FCS files for offline analysis in Kaluza, FCSEXPRESS, FlowJo, and other platforms.

STANDARDIZATION

Daily QC beads or any other reference material that is relevant for your application may be used as the standardization sample to set target values and calibrate the gain settings automatically.

LANGUAGE

English and Chinese

OPERATING SYSTEM

Windows® 7 Professional 64-bit

Windows® 8 Professional 64-bit

Windows® 10 Professional 64-bit

FCS FORMAT

FCS 3.0

MINIMUM COMPUTER SPECIFICATIONS

CPU: Intel Core i7, up to 3.9 GHz

RAM: 8 GB, 5 USB 2.0 and above ports

Storage: 1 TB drive in RAID1

Ethernet: Integrated 10M/100M/1000M GB

COMPENSATION

Automatic full matrix compensation

Manual full matrix compensation

Novel Compensation Library: store fluorescent spillover values of dyes to easily determine the correct compensation matrix with new gain settings

Import/export compensation values between experiments

Absolute linear gain amplification enables the use of compensation settings between experiments and sample types

INSTALLATION

DIMENSIONS (W X D X H)

Cytometer (with or without Plate Loader)

60.5 cm x 73.3 cm x 45.1 cm

Tanks and Holder

25 cm x 25 cm x 25 cm

WEIGHT

Cytometer: 79 kg

Cytometer with Plate Loader: 83.6 kg

POWER SPECIFICATIONS

Voltage: 100-240 V Power: 150 -250 W

OPERATING TEMPERATURE NON-CONDENSING

15-30 °C



For research use only. Not for use in diagnostic procedures.

© 2018 Beckman Coulter Life Sciences. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For Beckman Coulter's worldwide office locations and phone numbers, please visit "Contact Us" at beckman.com

FLOW-380006.18